


Course Overview	1
Unit 1: Introduction to Programming	9
Section 1 – What is LIMS Basic?	11
What is a variable?	14
What is a subroutine?	15
Variable types	16
What is a function?	20
Concatenating Variables	22
Section 3 – Some simple Functions	27
Date() and Time()	36
Chr()	34
InputBox()	31
Section 4 – Statements	39
What is a statement?	40
IF THEN	41
SELECT CASE	45
Loops	48
FOR NEXT	51
Section 2 – Using Functions and Statements in LIMS Programs	55
Status and Arguments	57
Section 1 – “LIMS” Functions	59
Commit Flag	60
Arrays in Function Arguments	61
Section 2 – Macro to Log samples	64
LogSample() and AssignTest()	65
Error Trapping	68
Progress Dialog	72
Macros	79
Section 3 – Running SQL queries from LIMS Basic	84
SQL	85

Table of Contents


Mr. R. L. Lee

The Query Window.....	68
The Stored Query Manager.....	91
The Data Explorer.....	93
The SQL() Function.....	95
Return arrays and Ubound()	98
Inserting variables in query string	104
Section 4 – Macro to close old projects	106
Section 5 – Macro to re-enable a user	112
Section 6 – Macro to write to a file	122
Introduction – Code triggers in LIMS	135
Section 1 – Event Triggers	137
What is an event trigger?	138
Context	141
Selecting from the context	143
Findings out what is in context	145
Breakpoint()	148
Calling a subroutine from an event trigger.....	150
Database event triggers	154
Section 2 – Automation Scripts	156
What is an automation script?	157
System variables	161
Adding a custom menu through an Automation script.....	162
Section 3 – Calculations	168
Section 4 – Template Formula	173
Section 5 – Format Calculations	178
Section 6 – ID Configurations	183
Final Note	186
Optional Exercise: Status Rule	187
Solutions to the Exercises	189

>>Duration: 10 minutes



- Domestic arrangements
- Restrooms



Objectives

LIMS experience

Personal background

Department / Functional
Area

Company

Introductions



1203-V6 - Overview - V02 © LabWare 2011 [3]

• Lab activity: QA, R&D, commercial?

• Personal: IT/database/programming experience?

• LIMS experience: experience or role

LIMS Basic Training:

Desired Learning Outcomes

1. Understand simple programming
2. Create subroutines using functions and statements within LabWare LIMS
3. Create Macros, Automation Scripts and Event Triggers using LabWare LIMS Basic functions
4. Recognise how the use of LabWare LIMS Basic may be useful in your functional area

1203-V6 - Overview - V02	[4]	© LabWare 2011
--------------------------	-----	----------------

Agenda

Agenda	
Day 1:	Introduction to LIMS Basic
Day 1:	• What is LIMS Basic, elementary concepts • Simple functions and statements
Day 2:	Using Functions and Statements in LIMS Basic
Day 2:	Running SQL from LIMS Basic Macros
Day 3:	Event Triggers, Automation Scripts, Calculations, Template Formulas, Format Calculations, Id Configuration
	1203-V6 - Overview - V02 © LabWare 2011
	[5]

1203-V6 - Overview - V02 [6] © LabWare 2011

- This is a "hands-on" course.
- Don't hesitate to explore and experiment with the interfaces as they are being explained to you. You won't break anything!
- If you are "stuck" please inform trainer for assistance.
- Ask any questions you like at any time.
- Please ensure cell phones are off or on vibrate.

Training Room

Course Overview

- Designed with both novice and experienced programmers in mind
- Pre-requisites:
 - LIMS Admin 1 course
 - A good all-around knowledge of LabWare LIMS configuration
 - LIMS Basic functions (developed and extended on an on-going basis)
 - Not all of the functions will be covered in this course
 - Additional resources / references:
 - On-line Help Files
 - LIMS Basic User Manual
 - LabWare White Paper: LIMS Basic Guidelines and Object Naming conventions (LabWare Support web page)

1203-V6 - Overview - V02


[7]

© LabWare 2011

- At the beginning of each exercise you will find a list of the functions and statements required to complete the exercise. At the beginning of each exercise you will find a list of the examples used in the exercises may not be the only way to achieve the same outcome. Solutions to the exercises are provided at the end of the training manual.

EXERCISES

1203-V6 - Overview - V02 [8] ©LabWare 2011



- This section is an introduction to LIMS Basic using simple variables, statements and functions.

By the end of this training segment you should be able to...

Objectives	
------------	---

1. Recognise and use different types of variables
 2. Create new subroutines
 3. Use simple LIMS Basic functions
 4. Understand the purpose of LIMS Basic statements
 5. Step through code and test the output

1203-V6 - Introduction to Programming - V02 © LabWare 2011

What is LIMS Basic?

Section 1

- LIMS Basic is a Scripting Language that is available in LabWare to allow users to introduce logic into the configuration of their system. For example:
- You could configure an Analyst's rule that uses LIMS Basic to automatically create an investigation and alert the QC Manager when the E. coli bug is found in a Microbiology sample.
- You could also configure an extra menu in the Project Manager that allows the user to re-allocate samples to another project.
- BASIC is the acronym for Beginner's All Purpose Symbolic Instruction Code. The BASIC programming language uses "ordinary English" to perform processor tasks. BASIC was introduced in 1964 and many variations of the language exist today.
- LIMS Basic is a variation of BASIC and is very simple to learn.
- LIMS Basic does not require variable typing, declaring or dimensioning as is the case with various other programming languages.
- LIMS Basic code is not compiled. The basic interpreter reads the code and translates the code in to actions. The code can be written and immediately tested.

1203-V6 - Introduction to Programming - V02
© LabWare 2011

4

LIMS Basic

- Scripting Language used by LabWare
- Variation of the BASIC programming language
- Simple scripting language to learn
- Interpreted Macro Language, not compiled

- All LIMS Basic code is stored in the LabWare Database (not in a separate executable)
- There are functions and statements that are similar to "BASIC" type statements, as well as LabWare LIMS specific functions
- Each statement must begin on a new line
- Statements are executed in the order that they appear (top to bottom)

1203-V6 - Introduction to Programming - V02 © LabWare 2011

5

Facts about LIMS Basic

- Code is stored in the LabWare database
- A mixture of "Basic" statements and functions, as well as "LabWare LIMS Basic" functions can be used
- Each statement begins on a new line
- Statements are executed in the order they appear

Elementary Concepts

Section 2

- Before the computer can perform the instructions or do the calculations it first has to be told exactly what to do. The instructions to the computer are called the PROGRAM.
- What we see, read or type for the program is simply for readability for ourselves, reviewers and users. The machine or computer only understands machine-readable code, such as binary (0's and 1's). Fortunately we do not have to do this translation ourselves; the programming language does this for us.
- We can easily define the answer as "7". The computer has to store each part of the sum within variables.
- A variable is assigned a value specifically or as a result of processing data.
- Variable values may be strings, numbers, arrays, date-times, times or dates.
- You may need to convert a variable's value type, e.g. from a number to a string to use the value in the MsgBox function.
- Variable names must start with a character. Variable names may not be a reserved word, e.g. USER, EMPTY, ERROR


- Instructions to a computer are called a PROGRAM
 - VARIABLES are placeholders used to store values
 - Example of a program to add two numbers:
- ```

Line 1 a = 2
Line 2 b = 5
Line 3 c = a + b


```
- The letters "a", "b" and "c" are VARIABLES
  - A variable can hold a value:
  - Specifically assigned to the variable
  - As a result of processing data
  - "Programming" allows a machine to perform tasks and calculations that people can do but in a more efficient manner.

## What is a variable?

- The easiest place to start writing or testing LIMS Basic code is in a **subroutine**.
- New Subroutines can be created in the Subroutines table.
- Subroutines are discrete pieces of code which can be called from various events within LabWare.



- Click on the **Test the Code** button on the toolbar. This will open the LabWare Debugger window.
- To execute the code line by line, click on the **Step** button.
- The first line of code above is not executed because it is a **comment**.
- Comments are added to code for clarification and documentation purposes.
- To “comment out” a line of code, simply make the line start with an apostrophe (').
- Comments can also be added at the end of a line.
- Comments can be added to a selection of code by highlighting the code and selecting the comment icon within the subroutine.



- The LabWare White Paper called "LIMS Basic Guidelines and Object Naming Conventions" provides best practices when writing LIMS Basic.
- While LIMS Basic is not case sensitive it is suggested that when identifying variables multi-case names are used for easier readability, e.g. sampleNumber or productName
- Tables and fields that are referenced in SQL calls or code are suggested to be in all CAPS
- It is important to add comments to the code, not only to describe what the code is meant to do, but to also describe the reason it was built to perform the way it was established.
- Meaningful comments are useful for reviewers, other programmers, and the support team to understand the purpose of the written code.
- Note: Not all lines of code need to be commented, just enough to add meaning.
- Code Snippets are defined block of LIMS Basic code (templates) that have been or can be created to assist with proper coding technique or to simplify coding.
- Code snippets are stored in the Code Snippets table.
- A block of code can be saved as a Code Snippet from the Formula Dialog window.
- A Code Snippet can be inserted into a Subroutine or Formula Dialog from the Formula window.

## Suggested Coding Syntax

- Use meaningful comments
- Use descriptive variable and array names
- Use indentation for nested code and loops
- Use UPPERCASE to identify TABLES and FIELDS
- Use code snippets when possible



## Suggested Coding Syntax

1203-V6 - Introduction to Programming - v02  
11 © LabWare 2011

In this exercise, you will create a subroutine with a simple calculation using variables.

1) Create a subroutine called NUMBERS.

2) Write a few lines of code which take the two values 7 and 9 and multiply them by each other.

3) Save the subroutine.

4) Test the LIMS Basic code by stepping through the code using the LIMS debugger. Notice the variable values in the right hand pane of the debugger window.

5) If time allows, add comments and code header to your subroutine.

**Exercise LB01-01: Variables**



- Although the “`diff = end - start`” code may look like an algebraic formula, it is not.
- It subtracts the variable values from each other and the answer is stored in the variable “`diff`”.
- Uses the two variables called “`start`” and “`end`” to store data or values.
- Example: `diff = end - start`
- As we saw in the previous section, variables can also be used to perform calculations. For example, “`diff = end - start`”.
- A variable does not have to be defined before it can be used
- Basic Manual.
- The variables can have any name, apart from a few reserved words listed in the LIMS Basic Manual.
- In LIMS Basic:
- The type of variable has to be specifically defined before it can be used.
- Variables have to have certain types of names depending on whether they are to store integers, floating point decimals, text, etc.
- The type of variable has to be specified before it can be used.
- In some programming languages:

1203-V6 - Introduction to Programming - V02  
© Labware 2011  
12

```
mass = 39.2932 = floating point (real number)
a = 2 = small integer
message = "Hello, this is a string" = text (string)
```

- Open the subroutine VAR\_TYPES to test this code.

- Here are some examples of values a variable might contain:

- LIMS Basic variables can have any name, apart from a few reserved words listed in the LIMS Basic manual.

## Variable Types

- However in programming languages, the value of  $x$  is redefined in each line as shown.
- Mathematical or algebraically speaking  $x$  would not be redefined in each line of this example. The value of  $x$  is always "0".
- 

1203-V6 - Introduction to Programming - V02      13      ©LabWare 2011


▪ Open the subroutine EQUATIONS to test this code.

|              |             |
|--------------|-------------|
| Value of $x$ | 0           |
| Line 1       | $x = 0$     |
| Line 2       | $x = x + 5$ |
| Line 3       | $x = x + 5$ |
| Line 4       | $x = x + 5$ |
|              | 15          |
|              | 10          |
|              | 5           |
|              | 5           |
|              | 10          |
|              | 15          |

▪ Although code may look like algebraic equations, programming languages do not follow algebraic behavior. Evaluate the following example:

**Variable Types**

- A LIMS Basic function is a pre-defined block of code that will perform a specific LabWare action.
- Think of all the things you can do interactively, e.g. Log a Sample, Add a Test, Enter a Result. There are literally hundreds of LIMS Basic functions. They are also known as commands.
- Most often LIMS Basic functions have arguments, which is the information contained in the parentheses following the function name. The values in the arguments can be used to perform the actions.
- Some LIMS Basic functions, for example arithmetic functions, are common to all BASIC languages. There are also many functions that are specific to the LabWare environment.
- TIP: It is strongly recommended that you always insert a LIMS Basic function. This will prevent syntax errors and reduce debugging.



- >>> Always insert, rather than type
- >>> Describe the use of the Quick Codes Selection Dialog

The screenshot shows a Windows application window titled "Quick Codes Selection Dialog". The window has a standard Windows title bar with icons for minimize, maximize, and close. The main area contains a list box with several items: "Unassigned", "Assigned", "Let", "FileFindReplace", "FormatWindowTextLeft", "RemoveSelSampleFromMatch", and "SampleValidValuesFromFiles". At the bottom of the window, there is a toolbar with three buttons: "Search Codes", "Full Search", and a magnifying glass icon labeled "Search Criteria". The status bar at the bottom displays the path "1203-V6 - Introduction to Programming - v02" and the number "15".

The magnifying glass icon in the toolbar is highlighted with a red box, and the text "Select Dialog" is written next to it.

- The Search Function icon can be used to open the Quick Codes Selection Dialog

**Inserting a LIMS Basic Function**

- In this example the Left LIMS Basic function is used: `part = Left(string,num)`
- `Left` is the name of the function.
- `part` is the return variable.
- `string` and `num` are the function's arguments; which are **separated** by a comma.
- The function processes the arguments and gives you an answer (value) in return.
- Create a new Subroutine called `TEST` to configure the code in the example.
- Test the code in the Subroutine.
- Notice the value of the variable `part` at the end of the code.
- What happens if you change the value of the variable `num` to 4?

- Example: The Left function can be used to extract a number of characters from a string
- Arguments are variables holding specific data values that will be used when the LIMS Basic function is executed
- The arguments can be defined as follows:

```
string = "LabWare"
num = 3
part = Left(string,num)
```


- The function will return the 3 left most characters of the string to the variable with the name `part`
- `part = "Lab"`

## Arguments in a LIMS Basic Function

- Message boxes can also be useful while building your code to detect errors.
- Message boxes are not only used to display a message to the user. As you will see later they can also be used to display a message to the user.
- The value of `string` is typed within double quotes to indicate that it is a string (`text`).
- `String` is argument in the function.
- The Message Box function: `MsgBox(string)`

1202-V6 - Introduction to Programming - V02 ©LabWare 2011

17



The screenshot shows a standard Windows-style message box. The title bar says "Information". The main text area contains the message "This text will appear on the screen" and features an exclamation mark icon. At the bottom right is a small "X" button, and at the bottom center is an "OK" button.

MsgBox(string)

`string = "This text will appear on the screen"`

- Example:
- `String` is the variable
- Syntax: `MsgBox(string)`
- Message boxes can be used to display "messages" to users.

**Function Demo: MsgBox()**

**Exercise LB01-02: Message Box**

In this exercise, you will create a subroutine with a simple string variable and you will display it in a message box.

None

Statements:

MsgBox(string)

Functions:


None

1) Create a subroutine called HELLO and insert the MsgBox function.  
2) Define the string variable as "Hello".  
3) Save the Subroutine.  
4) Test the code in the subroutine to display the variable's value.

1203-V6 - Introduction to Programming - V02  
© LabWare 2011  
18

- Through-out this training, you will notice that we have used either with no preference.
- It is a matter of preference to which you use.
- The & and + signs are inter-changeable to concatenate variables together.

1203-V6 - Introduction to Programming - V02  
19 ©LabWare 2011



MsgBox(string)

string = "The value of a is " + a + " and the value of b is " & b

b = 100

a = 50

Example:

The ampersand (&) or plus (+) sign is used to concatenate the values together.

Messages to be displayed can be built using a combination of strings and variables.

## Concatenating Variables

## Exercise LB01-03: Concatenating Variables



In this exercise, you will configure a subroutine that displays the first name and surname of John Smith in a message box.

**Functions:**  
None  
**Statements:**  
MsgBox(string)

1) Create a subroutine called USERNAME.

2) First create a string variable with the value "John".

3) Create another string variable with the value "Smith".

4) Create a third variable that concatenates both the first name and surname. Don't forget to insert a space between the two words.

5) Display the third variable in a message box.

20

© LabWare 2011

LB03-V6 - Introduction to Programming - V02

© LabWare 2011

## Exercise LB01-04: Concatenating more Variables



In this exercise you will configure a subroutine that calculates the circumference of a circle with a radius of 7 cm. You will display the result to the user in a message box.

Circumference formula:  $C = 2\pi r$ . Note:  $\pi = 3.14$

**Functions:**  
None  
MsgBox(string)

**Statements:**  
None

1) Create a subroutine called CIRCLE

2) Define each of the calculation parameters as separate variables.

3) Define a string variable that shows what the result is. Your string should be displayed like this:

4) Display the string in a message box.

The circumference of the circle is ..... cm

21

© LabWare 2011

1203-VB - Introduction to Programming - V02

1203-V6 - Introduction to Programming - V02  
© LabWare 2011  
22

## Some simple functions

### Section 3

1203\_V6 - Introduction to Programming - v02      23      © LabWare 2011

**Function Demo: InputBox()**

|                  |                                                                                                                                                                                                                      |
|------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Purpose:</b>  | The InputBox function displays a dialog box that prompts the user for input.                                                                                                                                         |
| <b>Syntax:</b>   | <code>x = InputBox(prompt, title, default)</code>                                                                                                                                                                    |
| <b>Comments:</b> | <code>prompt</code> is displayed in the dialog box.<br><code>title</code> is optional and is displayed in the title bar of the dialog box.<br><code>default</code> is default response if no other input is entered. |
| <b>Example:</b>  | <pre>prompt = "What is your name?"<br/>title = "NAME"<br/>default = ""<br/><br/>x = InputBox(prompt, title, default)</pre>                                                                                           |
| <b>Returns:</b>  | <code>x</code> contains the information typed by the user.                                                                                                                                                           |

- Use meaningful names for your variables.
- If you define the arguments inside the function you may save a few lines, BUT your code will be more difficult to debug.
- Spelling matters; copy and past or insert whenever possible! Double click on a word to highlight / select and then perform copy and paste.
- Do not try to type functions from memory. Use the Search Function icon to search for and insert a LIMS Basic function. This will help reduce error checking.

1. Insert functions with the Search Function icon to ensure syntax is correct.

2. Define variable arguments BEFORE the function and NOT inside it.

3. Give meaningful names to your variables.

|                                                                |                                                                                                                                                                        |                                                         |
|----------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------|
| <code>name = InputBox("What is your name?", "NAME", "")</code> | <code>x = "What is your name?"</code><br><code>y = "NAME"</code><br><code>title = ""</code><br><code>default = ""</code><br><code>prompt = "What is your name?"</code> | <code>n = InputBox(x,y,z)</code><br><code>z = ""</code> |
|----------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------|

**Good Programming Hints**

1203-V6 - Introduction to Programming - v02      25      ©LabWare 2011

Hello ..... so you are ... years old. That's interesting!

4) Display the result in a sentence like this:

3) Ask him to enter his age in a second input box.

2) First ask the user to type his name in an input box.

1) Create a subroutine called `GREETING_1`.

None

Statements:  
`MsgBox()`

Functions:  
`InputBox()`

In this exercise, you will write a subroutine that prompts the user to enter his first name and age in input boxes. The information will then be processed and displayed in a message box.

**Exercise LB01-05: Input Box**



- ASCII code definitions can be found online.
- The function for new line is Chr(10) and carriage return is Chr(13).
- Even the space bar and the carriage return have their own ASCII code. This allows them to be used in programming, for example a carriage return can be inserted so that it will display a string on a new line in a message box.
- ASCII code. For example, the ASCII code for the letter T is 0084. If you know the ASCII code of a character, you can usually insert it in a document by typing Alt+ASCII (Alt\_0233 will insert the letter e).
- Each character on your computer keyboard has its very own unique numeric code called an ASCII code. For example, the ASCII code for the letter T is 0084. If you know the ASCII code of a character, you can usually insert it in a document by typing Alt+ASCII (Alt\_0233 will insert the letter e).

1203-V6 - Introduction to Programming - V02      26      © LabWare 2011

**Chr()**

Purpose: The CHR function converts an ASCII code to its equivalent character.

Syntax: `x = CHR(integer)`

Comments: "integer" is an integer value in the range of 0 to 255. Refer to the table of ASCII characters for values.

Example 1:

```

 Returns: x = T
 x = CHR(65)

```

Example 2:

```

 Returns: Message = "This is the first line" + CHR(10)
 Message = Message + "This is the second line" + CHR(10)
 Message = Message + "This is the third line"

```

Returns: A message box displaying the following:

This is the first line  
This is the second line  
This is the third line

This is the third line

**Exercise LB01-06: Message Box**

In this exercise, you will write a subroutine that prompts the user to enter his first name and surname in input boxes. The information will then be processed and displayed in a multi-line message box.

Functions:

```

 None
 Statements:
 MsgBox()
 InputBox()
 Chr()

```

With multiple lines

The message box should look like this:

1) Create a subroutine called USER\_INFO.

2) First ask the user to enter his first name in an input box.

3) Ask the user to enter his surname.

4) Create a string variable to display the result in a message box on 2 separate lines.

| Function Demo: Date() |                                                    |
|-----------------------|----------------------------------------------------|
| Purpose:              | The Date function returns the current system date. |
| Syntax:               | today = Date()                                     |
| Comments:             | today is a date variable type.                     |
| Example:              | today = Date()                                     |
| Returns:              | today = 07/22/2000                                 |

## Function Demo: Time()

Purpose: The Time function returns the current system time.

Syntax: now = Time()

Comments: now is a time variable type.

Example: now = Time()

Returns: now = 11:09:45

1203-V6 - Introduction to Programming - v02      30      © LabWare 2011

**Exercise LB01-07: Date() and Time()**

In this exercise, you will create a subroutine with date and time variables and you will display them in a message box.

Functions:

- None
- Statements:

|          |        |
|----------|--------|
| Date()   | Chr()  |
| MsgBox() | Time() |

4) Notice the variable types as you test your code in the LIMS Debugger

3) Create a string variable that displays the date and time in a message box like this:

2) Insert the date() and time() functions.

1) Create a subroutine called DATE\_INFO.

Time: .....

Date: .....

## Statements

## Section 4

- LIMS Basic statements allow the introduction of logical structure within LIMS
- For example:

```

IF a certain condition is true (e.g. my sample is in spec)
THEN do one thing (e.g. auto-approve it)
ELSE do another thing (e.g. send an alert)
ENDIF

```

- These structures are very useful to automate "routine" decision making processes.
- Most LIMS Basic statements are generic BASIC statements, so if you have previous experience in this, you will find them quite easy to use.
- Once more, if you always insert statements instead of typing them in, you will save yourself a lot of syntax errors.

- LIMS Basic requires specific formats for some statements.
- For example, an IF THEN statement is a conditional structure. It allows you to introduce logic:

```

IF a certain condition is true (e.g. my sample is in spec)
THEN do one thing (e.g. auto-approve it)
ELSE do another thing (e.g. send an alert)
ENDIF

```

- Statements are lines of LIMS Basic code. Statements can consist of expressions, functions, and variables.
- LIMS Basic statements give structure to the code. You could say that if the functions were the "flesh" of the code, the statements would be the "bones".
- LIMS Basic requires specific formats for some statements.
- For example, an IF THEN statement is a conditional structure. It allows you to introduce logic:

1203-V6 - Introduction to Programming - V02  
©LabWare 2011  
32

• Available statements

```

IF a condition is true THEN
 do something,
ELSE
 do something else.
ENDIF

```

• Statements can be inserted using the Insert a Statement button [≡]

• Statements allow the introduction of logical structure within LIMS

• Statements allow the introduction of logical structure within LIMS

• Statements can consist of expressions, functions, and variables.

• Statements give structure to the code. You could say that if the functions were the "flesh" of the code, the statements would be the "bones".

• Statements are lines of LIMS Basic code. Statements can consist of expressions, functions, and variables.

• LIMS Basic requires specific formats for some statements.

• For example, an IF THEN statement is a conditional structure. It allows you to introduce logic:

```

IF a certain condition is true THEN
 do something,
ELSE
 do something else.
ENDIF

```

• Available statements

## What is a statement?

```

ENDIF
Msgbox(y)
ELSE
Msgbox(xy)
xy = x + , - , + y
y = Str(y)
x = Str(x)
ELSEIF (x = y) THEN
Msgbox(x)
IF (x < y) THEN

```

```

y = 10
x = 1

```

- For example:

- It is possible to test additional conditions using the ELSEIF.
- Optional ELSEIF**

1203-V6 – Introduction to Programming – V02      33      © LabWare 2011

**IF THEN Statement**

**Syntax:** IF (expression) THEN  
statements...  
ELSE  
statements...  
ENDIF

**Purpose:** The IF THEN statement allows conditional execution of statements based on the evaluation of a logical expression.


**Comments:** The expression must evaluate to True or False.  
The full expression must be enclosed in parentheses.

**Example:** x = 1  
y = 10

**Test this code in the IF THEN subroutine**

**ENDIF**  
MsgBox(y)  
ELSE  
MsgBox(xy)  
xy = x + , - , + y  
y = Str(y)  
x = Str(x)

- When testing or debugging LIMS Basic code the default behavior is not to step into a statement.
- Initially, the debugger highlights the whole statement. The step button will evaluate all conditions and execute the subsequent actions in one go.
- To step into the Statement click on the Step Into button. Once in the statement you can use either the Step or Step Into buttons to step through the code.
- Use the Step Out button to step out of the statement.



- Try using this in the next exercise instead of an `InputBox()`. What other functions might you need to determine if it is morning or afternoon?
- This is an additional function that can be used in place of prompting users with an `InputBox()`.

**Function Demo: `PromptForDateTime()`**

|                                                                                                                                                                                                                                                                                                                                                                                                    |                                                                                 |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------|
| <p><b>Purpose:</b> The <code>PromptForDateTime</code> function prompts the user for date and time value.</p> <p><b>Syntax:</b> <code>dt = PromptForDateTime(title, default)</code></p> <p><b>Comments:</b> <code>dt</code> is a <code>Date/Time</code> data type.</p> <p><b>Example:</b> <code>dt = PromptForDateTime()</code></p> <p><b>Returns:</b> <code>dt = 08/28/2010 01:24:36 PM</code></p> | <p>1203.V6 - Introduction to Programming - V02<br/>© LabWare 2011</p> <p>35</p> |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------|

- `PromptForDateTime()` will need the use of `Left()`, `Right()` and/or `Mid()`

1) Create a subroutine called `GREETING_2`

```

 IF THEN
 Statements:
 MsgBox()
 Functions:
 InputBox()
 MsgBox()

```

In this exercise you will work out whether it's morning or afternoon and you will display your finding in a message box.

**Exercise LB01-08: Finding out if it's morning or afternoon**



1203-V6 - Introduction to Programming - V02      ©LabWare 2011

36

2) First prompt the user to enter a time in a format HH:MM:SS (24 hours format)

3) Now work out whether the entered time is before or after 12:00:00

4) If it's morning, display the message: **Good morning!**

5) If it's afternoon, display the message: **Good afternoon!**

If time allows, try using the `PromptForDateTime()`

**SELECT CASE Statement**

Purpose: The SELECTCASE statement executes one of several statement blocks depending on the value of a variable.

Syntax:

```

SELECTCASE variable
CASE comparison list (e.g. "A", "B")
CASE comparison list (e.g. "C")
CASE statements...
CASE ELSE
statements...
ENDSELECT

```

Example: X = 3

Test this code in the subroutine

```

SELECTCASE x
CASE 2.4.6
MsgBox("x is even")
CASE 1.3.5
MsgBox("x is odd")
CASE ELSE
MsgBox("something else")
ENDSELECT

```

1203-V6 - Introduction to Programming - Vol2  
© LabWare 2011

37

**Right()**

**Purpose:** The Right function returns the rightmost number of characters in a string.

**Syntax:** `x = Right(string, num)`

**Comments:** If num is greater than the length of the string, the entire string is returned. If num is 0, an empty string is returned.

**Example:** `string = "LabWare"`  
`num = 4`  
`x = Right(string, num)`  
`x = "Ware"`

**Returns:** `x = "Ware"`

1203-V6 - Introduction to Programming - V02      38      © LabWare 2011

**Function Demo: Right()**

1203-V6 - Introduction to Programming - V02      39      ©LabWare 2011

4) If the number is even, display the message: <entered number> is even.

3) If the number is odd, display the message: <entered number> is odd.

2) Look if that number is odd or even. Tip: look at the last digit of the entered number.

1) Prompt the user for a number

**SELECTCASE**

Statements:

- Rig ht()
- MsgBox()
- InputBox()
- InputBox()

Functions:

In this exercise you will modify the SELECTCASE subroutine so that the user is prompted for a number to test. The subroutine should work for ANY number, no matter how big.

Exercise LB01-09: SELECTCASE



- Loops allow the repetition of a section of code for a number of times. For example, you might want to log several samples at once or test a series of conditions one after the other.

Loops execute a section of code a number of times

## LOOPS

There are several types of loops:

- DO LOOP
- WHILE WEND
- FOR NEXT

- condition is true

- The DO LOOP and WHILE WEND loops continue while a

- The FOR NEXT loop continues for a defined number of loops.

- `Sqr()` is the square root function

1203-V6 - Introduction to Programming - V02      © LabWare 2011

41

## DO LOOP Statement

**Purpose:** The DOLOOP statement executes a block of statements while a condition is true.

**Syntax:** DO WHILE (expression)

**Comments:** The expression (condition) logically evaluates to True or False. The statements are repeated while the expression is True. DO LOOPS may be nested.

**Example:** `x = 1`

**Test this code in the**

DO WHILE (`Sqr(x) < 5`)  
  `X = X + 1`  
  `MsgBox(X)`  
LOOP

**Subroutine:**

DO LOOP

**WHILE WEND Statement**

**Purpose:** The WHILE WEND statement executes a series of statements in a loop as long as a given condition is true.

**Syntax:** WHILE (expression) statements... WEND

**Comments:** The expression (condition) must be in parentheses.

The statements are repeated while the expression is True.

**Example:**    `x = 1  
y = 10  
WHILE (x < y)  
 msgBox(x)  
 x = x + 1  
WEND`

Test this code in the subroutine

- **TIP:** Make sure to step into ANY loop while testing using the LabWare Debugger
- WHILE WEND loop
  - Statements are executed until the condition is no longer valid
- DO LOOP
  - Statements are executed until the condition is no longer valid
- FOR NEXT loop
  - Statements are executed a defined number of times

1203-V6 - Introduction to Programming - V02      ©LabWare 2011

43

## FOR NEXT Statement

**Purpose:** The FOR NEXT statement executes a group of statements a specified number of times.

**Syntax:** `FOR i = 1 TO num STEP step`

**Comments:** The *i* is a numeric variable to be used as a counter.

The *num* is the final value of the counter.

The *STEP* is the increment for each loop. If not specified, it is assumed to be 1.

The statements between the FOR and NEXT are the program lines between the FOR and NEXT.

**Example:** `FOR i = 1 TO 10 STEP 2`

**Test this code in the subroutine**

`FOR i = 1 TO num STEP step`  
`statements...`  
`NEXT`

**Subroutine:**

`MsgBox(x)`  
`NEXT`

- FOR NEXT loops can be used to build multi-line messages.
- You will need to initialize your variable before the loop. In this case we do not want any extra spaces so we are just using double quotes with no space in-between.
- Every time the code runs through the loop, "i" is increased by one (STEP is not defined therefore default STEP 1).

1203\_V6 - Introduction to Programming - V02

© LabWare 2011

44

**FOR NEXT Example**

```

message = ""
FOR i = 1 to 10
 message = message + "This is line number " + i + Chr(10)
NEXT
MsgBox(message)

Test this code in the
FOR_NEXT_2
subroutine

```

**Exercise LB01-10: FOR NEXT**



In this exercise, you will write a subroutine that builds a "3" times table. Display the times table in a message box.

functions:

- Chr()
- MsgBox()
- Statements:
- FOR NEXT

like this:

1) Create a subroutine named TIMES\_TABLE

2) Use the code example given in the previous page to build the message box. The message box should look like this:

Information

OK


10 x 3 = 30  
9 x 3 = 27  
8 x 3 = 24  
7 x 3 = 21  
6 x 3 = 18  
5 x 3 = 15  
4 x 3 = 12  
3 x 3 = 9  
2 x 3 = 6  
1 x 3 = 3

etc.

2 x 3 = 6  
1 x 3 = 3

45 1203-V6 - Introduction to Programming - v02 © LabWare 2011

Notes ...



- New statements
  - RETURN
  - GOSUB
- New functions
  - FileWrite(), UpdateDialog(), LogSample(), Wait(), UpdateTableFields()
  - SQLSelect(), Str(), ExecuteQueryTag(), FileClose(), Ubound(), FileExists(), FileNew(), OpenDialog(), DateTime(), PromptForYesNo(), Save(), SQL(), DateTimAdd(), AssignTest(), Mid(), CloseDialog(), CloseProject(), DateAdd(), ODBCDateTimestamp()
- LIMS concepts
  - Comments and documenting code
  - Error trapping
  - Arrays
  - SQL queries
- This section introduces more LIMS Basic functions, including the SQL function.

1203-V6 - Using Functions and Statements - V02 [2] © LabWare 2011

By the end of this segment, you should be able to:

**Objectives**

1. Configure typical LIMS Basic functions
2. Execute SQL queries from a LIMS Basic subroutine
3. Understand arrays and use them in programming
4. Error trap your code to avoid unexpected errors
5. Comment and document your code

# “LIMS” Functions

## Section 1

[3]

©LabWare 2011

1203-V6 - Using Functions and Statements - V02

©LabWare 2011

- We have already seen some functions that are common in BASIC languages.
- Some LIMS Basic function categories also contain:
- Functions that are specific to LabWare objects such as a Batch, Folder, Lot, Sample, Test, and Result.
- These functions perform LIMS operations such as updating a batch, logging a sample, and assigning tests.
- Many of these commands involve making changes to the data in the database.

- 1203-V6 - Using Functions and Statements - V02 [4] © LabWare 2011
- LIMS Basic functions are organised in categories - e.g.
  - LIMS Basic functions are common to "BASIC" languages (Str, Left, Right, Std, Avg, Date, etc. ...)
  - Many of these functions are common to "BASIC" languages
  - Date functions
  - String functions
  - Arithmetic functions
  - Others are specific to LabWare, for example functions to log samples or assign tests.

## LIMS Basic Functions

- status**
- Status is a variable that returns either **True** if the function was successfully executed or **False** if it failed.
- arguments**
- LIMS-related functions often have numerous arguments; some optional and some mandatory.
- commitFlag**
- The commit flag argument indicates whether the changes are committed to the database or whether they are kept in memory to be saved later.
- Commit Flag = T**
- Changes are immediately saved to the database when the function is executed.
- Commit Flag = F**
- Changes are kept in memory to be saved later, i.e. at the:

  - Next system commit
  - Next save() function
  - When the next function with the commitFlag = T is executed
  - Keeping the changes in memory to be saved later is equivalent to a user clicking the cancel button at result entry.
  - This allows the LIMS Basic to be configured with checks and error traps that will allow for scenarios where specific changes should not be written to the database.

1203.V6 – Using Functions and Statements - V02 [5] © Labware 2011

```
status = function(arg1,arg2,arg3,...,commitFlag)
```

Most LIMS-related functions have the following structure:

**“LIMS” Functions**

|                                                                                                                                        |                                                                          |
|----------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------|
| status:                                                                                                                                | a variable returning <b>True</b> (= success) or <b>False</b> (= failure) |
| function:                                                                                                                              | the function's name                                                      |
| arg(n):                                                                                                                                | the function's arguments (or inputs)                                     |
| commitFlag: a boolean flag ( <b>T</b> or <b>F</b> ) to indicate if the changes must be immediately committed (= saved) to the database |                                                                          |

- Functions that do not have a commitFlag:
  - The only way to find out is to reference the LIMS Basic Manual
  - May either auto-commit or not commit
- Functions that are not committed (by choice or not) must be committed later in the code with the Save() function or one of the other ways listed on the previous page.

1203-V6 - Using Functions and Streamers - V02 [6]

**To Commit or Not Commit**

- Some reasons why you may want to use commitFlag = "F":
  - Be able to opt out of committing the data
  - Commit all changes from a loop, in one operation after the loop
  - The Save() function will commit all uncommitted data.
- Some older functions in the "LIMS" category do not have a commitFlag
  - Others won't (e.g. AssignTest) therefore must be committed separately with the Save() function.
  - Some of these functions will auto-commit (e.g. LogSample)

- Refer to the LIMS Basic Manual to determine whether the array type argument is a one dimensional or multi dimensional array.

#### Important Tip:

- LIMS Basic code can be used to assign values to the elements.
- An array can consist of one or many elements.
- Two types of arrays:
  - One dimensional array (e.g. `arrayName[1] = "a"`)
  - Multi dimensional array (e.g. `arrayName[1,1] = "a"`)
- Some LIMS Basic functions may have **array** type arguments.
- An array is a variable that may contain more than one value, i.e. a data set

## Arrays

- The different values in the array are called array elements. In LIMS Basic an array can have as many elements as necessary.
- The array elements are defined using an index number within square brackets after the array name.
- Values can be assigned to the array elements

- 1203-V6 - Using Functions and Subroutines - V02      ©LabWare 2011
- Example of a one dimensional array:
  - Think of a one dimensional array as a list of values.
  - Elements of a single dimensional array are defined by a single index number, e.g. [1].
  - Elements of a multi dimensional array are defined by using one index number, e.g. [1,1].
  - Elements of a multidimensional array are each dimension, e.g. [1,1,1].
- 
- ```

graph TD
    A[Name of the array] --> B[Index for the array element]
    B --> C[Value assigned to the array element]
    
```
- ```

arraySAMPLENumbers[2] = 87543
arraySAMPLENumbers[1] = 87543
arraySAMPLENumbers[0] = 87544
arraySAMPLENumbers[3] = 87545

```

## One Dimensional Arrays

## One Dimensional Array: Sum()

- Function:  $x = \text{Sum}(\text{array})$
- Returns the sum of the values in the array
- Define the array elements and assign values
- Think of the one dimensional array as a list of array elements
- Each element can hold a value

|    |
|----|
| 10 |
| 29 |
| 39 |

```
array[1] = 10
array[2] = 29
array[3] = 39

x = Sum(array)
```

## Macro to log samples

## Section 2

[10]

1203-V6 - Using Functions and Statements - V02

© LabWare 2011

- The `LogSample()` function includes an automatic commit.
- `valuesArray`; as long as values are valid for the selected template.
- There is no limit on the number of array elements that can be defined in the `fieldsArray` and `valuesArray`.
- `valuesArray` is a one dimensional array that holds the values that will be used to populate the corresponding fields on the Sample Login Template.
- `fieldsArray` is a one dimensional array that holds the Sample Login Template fields that will be populated with values.
- `templateName` is the variable that holds the name of the Sample Login Template to be used to log the Sample.
- `sampleNumber` is the variable that will hold the new Sample Number.
- The `sampleNumber = LogSample(templateName, fieldsArray, valuesArray)` function logs a Sample and returns the new Sample Number.

1203-V6 - Using Functions and Structures - V02 [1] © LabWare 2011

```

sampleNumber = LogSample(templateName, fieldsArray, valuesArray)
valuesArray[2] = "CHOC-CAKE"
valuesArray[1] = "ALPHA"
fieldsArray[2] = "PRODUCT"
fieldsArray[1] = "CUSTOMER"
templateName = "CUST-PROD"

```

• `valuesArray` lists the values to put in the fields when logging the sample  
 • `fieldsArray` lists fields in the "CUST-PROD" Sample Login Template

• `sampleNumber = LogSample(templateName, fieldsArray, valuesArray)`  
 • Logs a Sample and returns the new Sample Number (auto commit)  
 • `sampleNumber = LogSample(templateName, fieldsArray, valuesArray)`

**Function Demo: LogSample()**

- The AssignTest function assigns a Test to a Sample.
- SampleNum - the Sample Number of the Sample (Test) to be assigned to the Sample
- analysisName - the name of the Analysis (Test) to be assigned to the Sample
- numReps - the number of Test replicates to be assigned to the Sample
- fieldsArray - a one dimensional array containing the Test record field names to be updated by the function
- valuesArray - a one dimensional array with values with which the new Test records will be assigned to the Sample
- maxReps - the maximum number of Test replicates to be assigned to the sample
- parentTestNumber - Test number of the Test to act as a Parent Test, if null, the Test is assigned to the Sample
- Note also that there is no database commit as part of the AssignTest() function. You will need to use the Save() function to save the new tests to the database.

1203-V6 - Using Functions and Statements - V02 [12] © LabWare 2011

```

Example:
parentTestNumber,
"testNumberArray", fieldsArray, valuesArray, maxReps,
status = AssignTest(sampNum, analysisName, numReps, fieldsArray[1] = "TEST LOCATION",
 fieldsArray[2] = "TEST COMMENT",
 fieldsArray[3] = "PH",
 analysisName = "PH",
 numReps = 3,
 sampNum = 741
 valuesArray[1] = "BGM",
 valuesArray[2] = "Test requested by customer",
 valuesArray[3] = "TEST COMMENT"
 status = AssignTest(sampNum, analysisName, numReps,
 "testNumberArray", fieldsArray, valuesArray, maxReps,
 fieldsArray[1] = "TEST LOCATION",
 fieldsArray[2] = "TEST COMMENT",
 fieldsArray[3] = "PH",
 analysisName = "PH",
 numReps = 3,
 sampNum = 741
 valuesArray[1] = "BGM",
 valuesArray[2] = "Test requested by customer",
 valuesArray[3] = "TEST COMMENT"
 status = AssignTest(sampNum, analysisName, numReps,

```

- This function can be used to assign a Test to a Sample

## Function Demo: AssignTest()

## Exercise LB02-01: Log a sample and assign tests

In this exercise, you will configure a subroutine to log a sample and auto-assign a test.



| Functions:                                                                                                                                        |             |
|---------------------------------------------------------------------------------------------------------------------------------------------------|-------------|
| AssignTest()                                                                                                                                      | LogSample() |
| AssignSample()                                                                                                                                    | Save()      |
| None                                                                                                                                              | Statements: |
|                                                                                                                                                   |             |
| 1) Create a subroutine called MACRO_LOG_SAMPLE.                                                                                                   |             |
| 2) Use the LogSample() function to log a single sample using the CUST_PROD values:                                                                |             |
| 3) Use the AssignTest() function to auto-assign the MOISTURE analysis in duplicate to the sample. You will also update the following TEST fields: |             |
| 4) Now test your code carefully.                                                                                                                  |             |

1203-V6 - Using Functions and Statements - V02 [13] © LabWare 2011

- Many of the LIMS Basic functions return a value of either true or false to the **status** variable.
- Furthermore, if the returned status is "False", LabWare provides a system variable called **lastError** that contains a high level indication / explanation of the error encountered.
- The status and lastError values can be used for error trapping.
- The value of the lastError variable:

  - Can be displayed to the user at the time that the error is encountered or it can be recorded in a log.
  - Will be displayed in the LIMS debugger window if you are using it to test your code.
  - Will be reflected in the debug.log file.

1203-V6 - Using Functions and Statements - V02      (14)      ©LabWare 2011

**Error Trapping**

Most functions can be error trapped based on the value returned for the **status** variable.

**status = function(arg1,arg2,...)**

- **status = true:** the function was executed successfully
- **status = false:** the function has failed
- **lastError:** system variable defined by LabWare when a LIMS Basic function failed

**Exercise LB02-02: Error trapping**



In this exercise, you will update the MACRO LOG-SAMP subroutine to include an error trap on the test assignment.

Additional functions:

Statements:

```

IF THEN
 MsgBox()

```

1) Check for the return status of the AssignTest() function.

If the status is false, display a message to the user that indicates what the error is. You may use the lastError system variable. Your message box should read:

If the tests have been successfully assigned, the message box should read:

Tests successfully assigned to sample .....  
2) Test your code.

1203-V6 - Using Functions and Statements - V02 [15] © LabWare 2011

**Exercise LB02-03: Log several samples in a loop**



In this exercise, you will create a new subroutine to log several samples instead of just one. The code is similar to the previous exercise, but this time, you will start by prompting the user for a number of samples to log.

Additional functions:

```
FOR NEXT
InputBox()
```

Additional statements:

```
FOR NEXT
InputBox()
```

1) Save your previous subroutine as MACRO\_LOG\_SAMPS


2) Prompt the user for a number of samples to log using an Input Box.

3) Use this information to set up a FOR NEXT loop that will log each sample and assign the MOISTURE test in duplicate.

1203-V6 - Using Functions and Statements - V02  
© LabWare 2011 [16]

## Exercise LB02-04: Build a multi-line message box

In this exercise, you will update the MACRO LOG\_SAMPS subroutine to display a message with information on ALL the samples logged. It should look like this:



Additional statements:  
None

Additional functions:  
Chr()

Tests successfully assigned to sample 15  
Tests successfully assigned to sample 16  
Tests successfully assigned to sample 17

- 1) Remove all message boxes from the FOR NEXT loop.
- 2) Initialize the message BEFORE the loop (message = "")
- 3) For each logged sample, build one line of the message using the technique:  
message = message & "Tests successfully assigned to sample ...."
- 4) Add a message box AFTER the loop to display your multi-line message.

1203.V6 - Using Functions and Statements - V02 [17] © LabWare 2011

- Progress dialogs are very useful for sections of code that are included in loops.
- Progress dialogs allow users to monitor the progress when the code is running thus users are less likely to believe that the system is "hanging".
- Progress dialogs consist of three functions
  - `OpenProgressDialog()` initializes the dialog
  - `UpdateProgressDialog()` is placed inside the loop and updates the progress
  - `CloseProgressDialog()` closes the dialog

1203-V6 - Using Functions and Subroutines - V02 [18] © LabWare 2011

**PROGRESS\_DIALOG**

You can test the above dialog in the subroutine

Square root of 4 is 2.0  
Square Roots

- 3 functions are required:
- `OpenProgressDialog`
- `UpdateProgressDialog`
- `CloseProgressDialog`

Progress dialogs display the progress of an operation to the user and are often used in loops.

## Progress Dialogs

- allowCancelFlag** – a boolean flag that indicates whether the user is allowed to cancel the current operation
  - showPercentage** – if “T” (true) a % value and progress bar will be displayed in the dialog window
  - message** – the displayed message in the dialog window
  - title** – the displayed name in the title bar of the dialog window
- OpenProgressDialog()**

1203-V6 – Using Functions and Statements - V02 [19] © LabWare 2011

**Comments:** The allowCancelFlag is an optional argument.

**Syntax:** status = OpenProgressDialog(title,message,showPercentage,allowCancelFlag)

**Purpose:** The OpenProgressDialog function displays a progress dialog.

**Example:**

```
status = OpenProgressDialog("Progress Dialog", "Testing 1, 2, 3", "T", F)
```

**Returns:** status = true if the function was successful, false otherwise.

**allowCancelFlag** = “F”  
**showPercentage** = “T”  
**message** = “Testing 1, 2, 3”  
**title** = “Progress Dialog”

**OpenProgressDialog()**

**Function Demo: OpenProgressDialog()**

- percentage displays a percentage value and controls the progress bar indicator display
- message is displayed in the dialog box

## UpdateProgressDialog()

1203-V6 - Using Functions and Subroutines - V02 [20] ©LabWare 2011

**Purpose:** The UpdateProgressDialog function updates the message and progress bar in a progress dialog.

**Syntax:** status = UpdateProgressDialog(message,percentage)

**Example:**

```
status = UpdateProgressDialog("Testing 4, 5, 6", 20)
percentage = 20
message = "Testing 4, 5, 6"
```

**Returns:** status = true if the function was successful, false otherwise.

**Function Demo: UpdateProgressDialog()**

- It is good practice to add all three functions to your code at the same time.
- The progress dialog will not close automatically leaving an open function.
- It is very important to make sure that you close your progress dialog.

1203 V6 - Using Functions and Statements - V02 [21] © LabWare 2011

**CloseP**rogressDialo**g**

Have You  
Closed Your Progress Dialog?

**Syntax:**    status = CloseProgressDialog()

**Purpose:**    The CloseProgressDialog function closes a progress dialog.

**Returns:**    status = true if the function was successful, false otherwise.

**Function Demo:** `CloseProgressDialog()`

1203-V6 - Using Functions and Statements - V02 [22] © LabWware 2011

## Function Demo: Wait()

**Purpose:** The Wait function causes the system to pause for a specified period of time.

**Wait()**

**Syntax:** status = Wait(seconds)

**Example:**

```
status = Wait(3)
seconds = 3
```

**Returns:** status = true if the function was successful, false otherwise

1203-V6 - Using Functions and Statements - V02 [23] © LabWare 2011

Good Programming Reminders



- Comment your code
- Good comments make it easier to understand and debug code
- Use correct indentation
- especially within loops

In this exercise, you will update the MACRO LOG\_SAMPS subroutine, to include a progress dialog showing the progress of the sample logging.

**Additional Functions:** OpenProgressDialog(), UpdateProgressDialog(), CloseProgressDialog(), WaitForNone  
**Additional Subroutines:** None

1) Insert a line just before the FOR NEXT statement to open and initialize the progress dialog.

2) The updating of the dialog should take place inside the FOR NEXT loop, towards the end. Ensure that the percentage calculation is correct. You may want to add the Wait() function.

3) Don't forget to close the progress dialog at the end of the subroutine.

4) Add some comments. Make sure that the code inside the loop is correctly indented.

This subroutine has 3 distinct parts: prompting for a number of samples and initialize progress dialog, loop to log samples and assign tests, and close progress dialog and display message. Look at the exercise solution for tips on how to present your code.

## Exercise LB02-05: Progress Dialog



- Macros can be used to execute repetitive or frequently used tasks such as displaying reports or re-printing labels.
- Macros may also lessen the need to train users on LIMS functionality, for example if you create a sample login macro the user may not need to know how to log a sample, which template to select and whether to assign it to a project, lot or batch.
- Macros may also provide a way to provide custom workflows. Prior to visual workflows, macros were the only way to provide custom workflows.
- In LabWare v6 it is recommended that visual workflows be used to provide custom workflows.
- **Important:** Prior to visual workflows, macros were the only way to provide custom workflows.

Subroutine function

• Macros can call a Subroutine with the GoSub statement or the Make a Macro available

• The LabWare toolbar needs to be configured for a user's role to Macro records are stored in the User Programs table in LabWare.


• Macros are user-added menus that are available from the LabWare main menu.

**Macros**

The screenshot shows the LabWare software window. The title bar reads "123-V6 - Using Functions and Subroutines - V02" and "© LabWare 2011". The menu bar includes "File", "Configure", "Run", "Updates", "Debug", "Search", "Help", "Workflows", and "Macros". A dropdown menu for "Macros" is open, showing options: "(1) User dialog", "(2) User window", "(3) User report", and "(4) Lab grid". Below the menu bar, there is a toolbar with several icons. The main workspace is currently empty.


- Select the Explorer Table Manager icon on the LabWare Main Menu
- Select the USER PROGRAMS table
- Create User Program record with the name, "HELLO"
- Define the Menu title as "Hello" (This is the menu that will display in under the Macro Menu)
- Configure Your MsgBox("Hello")

## Creating a User Program



- Log out of LabWare, log back in
- Select Configure > Menus and toolbar...
- Assign a bitmap to the macro (if a bitmap will be used)
- Drag the new menu under Macros in the All Menus pane
- Drag bitmap onto Toolbar pane and click OK (if a bitmap will be used)
- If Macro is not displayed, log out and log back in

### Addling the Macros Menu item to the LabWare Main Menu



>> If time, demo how a visual workflow would be used for the same subroutine.

- We will discuss another function to call a subroutine, `Subroutine()`, in the next unit.
- Subroutines can be called from Macros using the `GoSub` statement.
- Using subroutines instead of placing code directly within your user programs allow you to create generic subroutines that may be used in multiple places.
- It also allows for more efficient change control. If your subroutine needs updating you only need to update one area instead of each individual user program.

1203-V6 - Using Functions and Statements - V02 [28] © LabWare 2011

**Example:** `GoSub MYSUB`

**Syntax:** `GoSub SubName`

**Purpose:** The `GoSub` statement is for execution of subroutines.

**Comment:** Use the "Insert Statement" icon to insert the `GoSub` statement.

- 1) Create a new User Program called LOG\_SAMPS. The menu title for the macro should be displayed as "Log samples".  
 2) Insert the GOSUB statement to call the MACRO LOG\_SAMPS subroutine and save.  
 3) Log out of LIMS and back in.  
 4) Configure the LabWare Main Menu to display the Macros menu item.  
 5) Test the Macro.

You can now configure a macro to execute the code that you have configured in the previous exercise.

### Exercise LB02-06: Set up a macro



LB02-06 - Using Functions and Statements - V02 [29] ©LabWare 2011


| Functions:  | None  |
|-------------|-------|
| Statements: | GOSUB |

# LIMS Running SQL from Section 3

**Reminder:** SQL stands for Structured Query Language. It is a generic language used to query or modify databases.

- LIMS Admin 1 training and that you may have used to build LabWare query tags.

- The following few slides provide a short review of SQL principals that were discussed in LIMS Admin 1 training and that you may have used to build LabWare query tags.
- For instance, you may want to list all the grades for a product in a dropdown list, or you can add them to a project.
- There are many instances in LIMS Basic when you will need to query the database and process the result of that query.
- For instance, you may want to find out which samples belong to a certain location so that you may want your code to find out which samples belong to a certain location so that you can add them to a project.



- Data regarding customers is stored as individual records in the table, e.g. ALPHA, BETA, etc.
- The name of the table is CUSTOMER
- The above table is an example of a database table.

1203-V6 - Using Functions and Subqueries - V02 [32] © LabWare 2011

The screenshot shows the Microsoft Access application window titled "Customer". The main area displays a table with the following data:

| ID    | NAME  | GROUP_NAF      | DESCRIPTION                                     | COMPANY_I      | ADDRESS    | ZIP   |
|-------|-------|----------------|-------------------------------------------------|----------------|------------|-------|
| ALPHA | ALPHA | DEFAUL         | The Alpha com Alpha Ltd                         | Aspen Road     | Aberdeen   | UK    |
| BETA  | BETA  | DEFULT         | The Beta comp Beta ltd                          | Double Road    | Bromleyham | UK    |
| SOFA  | SOFA  | PPC CC DEFUALT | The Soda Pop Co The Soda Pop Co The Soda Pop Co | 123 Bubbles Av | Gastond    | CA 12 |

The status bar at the bottom of the screen displays "Customer".

- This slide is an example of how a SQL query works. The selected data (or query result) will be returned in a multi-dimensional array.
- Using **Select \* from (TABLE)** returns all data from that table. This is useful when you are not sure what kind of data is stored in each of the table's fields. For example, to determine how the ADDRESS\_3 field is being used.
- When selecting a FIELD from a TABLE, the return will be an array.
- If you specify a WHERE clause, the text value must be in single quotes.

**Selecting Data from a Table**

| CUSTOMER table |               |            |           |       |
|----------------|---------------|------------|-----------|-------|
| NAME           | ADDRESS_1     | ADDRESS_2  | ADDRESS_3 |       |
| ALPHA          | Aspen Way     | London     | England   | GAMMA |
| BETA           | Boston Avenue | Birmingham | England   | BETA  |
| Gamma          | Greenway Road | Glasgow    | Scotland  | ALPHA |

Examples of SELECT queries:

```

Select * from CUSTOMER;
Select NAME from CUSTOMER;
Select NAME from CUSTOMER where NAME = 'BETA';
Select ADDRESS_2 from CUSTOMER where NAME = 'BETA';
Select NAME, ADDRESS_2 from CUSTOMER where ADDRESS_3 =
 'Birmingham';
Select * from CUSTOMER where ADDRESS_3 = 'England';

```

[33] 1203-V6 - Using Functions and Statements - V02 ©LabWare 2011

Where FIELD<sub>4</sub> = numeric value and FIELD<sub>5</sub> = 'text value'  
 From TABLE\_NAME  
 Select FIELD<sub>1</sub>, FIELD<sub>2</sub>, FIELD<sub>3</sub>....


- Select SQL queries generally follow the syntax:

1203-V6 - Using Functions and Statements - V02 [34] © LabWare 2011

**SELECT Query**

SELECT NAME, ADDRESS<sub>2</sub>  
 FROM CUSTOMER  
 WHERE ADDRESS<sub>3</sub> = 'England'  
 -The table  
 -The selected fields  
 -The criteria that records  
 have to match to be selected  
 -Text values in the WHERE clause must be framed with single quotes

- Important: Only execute SELECT type SQL Queries from the LabWare Query Window
- The Structure button will display tables and fields. In the next few slides we will look at other areas queries can be written.
- The Query Window is a useful tool when developing LIMS Basic.



The screenshot shows the LabWare Query Window titled "1203-V6 - Using Functions and Statements - V02". The window has a menu bar with "File", "Edit", "View", "Format", "Tools", "Help", and "Query". A toolbar below the menu includes "New", "Open", "Save", "Print", "Exit", "Copy", "Paste", "Find", "Replace", "Select All", and "Clear". The main area contains a table with columns "ID", "NAME", "DESCRIPTION", and "SELECT FROM LOCATION". The table lists several rows of data. At the bottom of the window, there is a status bar with the text "Query window".

Demo: Query Window

- SQL queries can be executed from the Query Window
- Select Debug > Query Window to open the window.
- Type query in the top pane.
- The selected data is returned in the bottom pane by clicking on the Showit button.

## Exercise LB02-07: The Query Window

In this exercise, you will create and run a SQL query from the Query Window.

- 1) Open the Query Window and configure a query that will:
  - select the SAMPLE NUMBER
  - TEXT\_ID
  - PRODUCT
  - where the PRODUCT value was defined as CHOC\_CAKE
  - from the SAMPLE table for samples
  - Click on the ShowAll button to display the data.

[36]

LB02-06 - Using Functions and Subqueries - V02

LB02-01 - Advanced Configuration Using LIMS Basic

- To execute a Stored Query click on Perform Query or select File > Run Query. Once the query is executed, the results can be exported to Excel.
- 1) File > Open Stored Query Manager
  - 2) Select New Query; Enter name for query
  - 3) Highlight the table you are selecting from. Click Next
  - 4) Click Add button to define fields you want to report on. Click Next
  - 5) Click Add button to define fields for the "where" clause. Click Next
  - 6) Define the values for the "where" clause. Click Next
  - 7) Define a sort order (optional)
  - 8) Limit number of returned records (optional)
  - 9) The resulting query is displayed. Click Finish

### Using Stored Query Manager

1203-V6 - Using Functions and Subqueries - V02 [37] ©LabWware 2011

1) Select File > Open Stored Query Manager

2) Queries can be created and stored with the Stored Query Manager Wizard

Demo: **Stored Query Manager**

- In this exercise, you will create and run a SQL query from the Stored Query Manager.
- 1) Open the Stored Query Manager and build the same query as in the previous exercise – i.e. select the SAMPLE NUMBER from the SAMPLE table for samples where the PRODUCT value was defined as CHOC CAKE
- 2) Save the query.
- 3) Click on the "Perform Query" button to display the data.

### Exercise LB02-08: The Stored Query Manager



1203-V6 - Using Functions and Statements - V02 [36] © LabWare 2011

1203-V6 - Using Functions and Statements - V02 [36] © LabWare 2011

- The Data Explorer allows users to export the results either to Excel or to an adhoc folder.
- Data Explorer allows the use of prompts for your „where“ clause.
- Select the fields first by table, then select the field from the drop down in each column/row
- Enter your „where“ clause
- Click on the Run icon to run the query
- To view the SQL select Query > Edit SQL
- Data Explorer and Store Query Manager are useful tools to assist with trending and quick information gathering that end users may be granted restricted access to.

### Using Data Explorer

1203-V6 - Using Functions and Subelements - V02      © LabWare 2011

The screenshot shows the Data Explorer window with the title "1203-V6 - Using Functions and Subelements - V02". The main area displays a list of records with columns for "Text Id" and "Text". Below the list is a toolbar with various icons for filtering, sorting, and manipulating data. The status bar at the bottom shows "Data Explorer Projects".

- File > Open Data Explorer
- The Data Explorer is another place in LIMS where queries can be created, stored and run from.

**Demo: Data Explorer**

1203-V6 - Using Functions and Statements - V02 [40] © LabWare 2011

In this exercise, you will create and run a SQL query from the Data Explorer.

1). Build the same query as in the previous exercise in the Data Explorer.

2) Save the query.

3) Click on the "Run Query" button to display the data.

**Exercise LB02-09: The Data Explorer**

